


MC7741 Emergency Shut-off Valve Box DN15. Product description.



Product Description

For safety and service reasons a central gas system must be equipped with shut-off valves placed so that the gas supply can easily be interrupted. The valves are mounted in a box. The emergency shut-off valve boxes shall be placed so that the gas can be shut off section wise. This means that the boxes should be placed before each ward, operating unit, part of ward for critical treatment and individual surgeries.

The emergency shut-off valve box is delivered with connection tubes and each box has been test pressurized and controlled for tightness. The emergency shut-off valve has large ergonomical handles.

If mounted in a recessed way, the emergency shut-off valve box fits walls with 70 mm beam. With a 90 mm beam there is extra space (23,5 mm) behind the valve box usable for e.g. fire isolation.

All models, also with four or five gases, fit between the beams in a CC-60 wall. The box is gas-tight which prevents gas accumulation inside the wall.

It is important that the boxes are placed so that they are easily available for authorized personnel. The front glass shall be sealed.

In order to avoid mistakes the boxes shall be clearly and distinctly marked with gas sort. A sign showing which section the box serves must be placed in its immediate vicinity.

The valves are open when the handles are in vertical position in a line with the printed marking on the plate.

Recycling

The product shall be recycled according to local regulations.

Technical Datas

Maximum working pressure	16 bar
Tube dimension	diameter 15 mm
Pressure class	PN 16
Regulatory status	Complies with medical devices directive 93/42/EEC. Complies with EN ISO 7396-1 and present SIS HB 370.
Classification	Class IIa
CE-marking	CE 0434

Article Numbers**Emergency shut-off valve box with turning copper tubes**

Art. No.	Denomination
321561	1 pc valve DN 15
321562	2 pcs valve DN 15
321563	3 pcs valve DN 15
321564	4 pcs valve DN 15
321570	5 pcs valve DN 15

Accessories

Art. No.	Denomination
321583	Frame for outside mounting for box for 1-3 valves, turning
321584	Frame for outside mounting for box for 4-5 valves, turning

Emergency shut-off valve box with through copper tubes

Art. No.	Denomination
321571	1 pc valve DN 15
321572	2 pcs valve DN 15
321580	3 pcs valve DN 15
321581	4 pcs valve DN 15
321582	5 pcs valve DN 15

Accessories

Art. No.	Denomination
321585	Frame for outside mounting for box for 1-3 valves, through
321586	Frame for outside mounting for box for 4-5 valves, through

Contact us

Linde Healthcare is committed to working with healthcare providers and regulatory authorities to continuously promote safe use of medical products and improve patient care.

We provide medical gas products, therapies, technical solutions and services to hospitals, clinics, nursing facilities, emergency management services and home healthcare providers around the world.

With our long experience and understanding of healthcare realities, you can depend on solutions that are delivered and serviced to the highest possible standards of quality, safety and efficiency.